

EpiCypher[®]

Bringing Epigenetics to Life



Strattech

CUTANA[™]

Ultra-sensitive genomic
mapping assays for
ChIC / CUT&RUN


CUTANA™ Products and Services

CUTANA™ is an ultra-sensitive genomic mapping assay that uses a proprietary immunotethering approach similar to **Chromatin ImmunoCleavage (ChIC)** and **Cleavage Under Targets and Release Using Nuclease (CUT&RUN)** methods.

CUTANA assays are used to map histone post-translational modifications (PTMs) and chromatin interactors with high resolution. This approach leverages a factor specific antibody to tether a fusion of protein A, protein G and micrococcal nuclease (pAG-MNase) to genomic binding sites in intact cells, which is then activated by the addition of calcium to cleave DNA.

FIGURE 1

Overview of the CUTANA approach. Cells (or nuclei) are immobilized on lectin-coated magnetic beads, permeabilized, and incubated with an antibody to a chromatin target (e.g. histone PTM or chromatin / DNA binding protein). Next, pAG-MNase is added and activated via Ca²⁺. The clipped chromatin fragments diffuse out, followed by DNA purification and next-generation sequencing. CUTANA assays vastly outperform those using chromatin immunoprecipitation, the current gold-standard genomic mapping assay.



CUTANA vastly outperforms leading ChIP-seq approaches

TABLE 1


Comparison of CUTANA assays with ChIP-seq. Low sequencing requirements, improved Signal : Noise, and reduced assay time make CUTANA far superior to ChIP-seq.

FEATURES	CUTANA™	ChIP-Seq
Cells required	100	10,000*
Protocol time (cells - DNA)	< 1 day	> 2 days
Sequencing reads required	3 million	30 million
Signal : Noise	High	Low

*Source: Dirks et al., Clin Epigenetics. 2016 Nov 21;8:122

CUTANA is proprietary technology based on ChIC (US20070009937A1; Schmid et. al, Mol Cell 2004) and CUT&RUN (patent pending; Skene and Henikoff, eLIFE 2017) methodology.

The CUTANA workflow - From cells to data in less than 4 days



CUTANA saves you nearly 5-10 fold in sequencing requirements.

CUTANA assays generate superior quality data compared to ChIP-seq using as little as 3 million sequencing reads.




FIGURE 2

A representative region on chromosome 20 of H3K27me3 from K-562 cells is shown, comparing standard ChIP-seq methods and sequencing depth (black) vs CUTANA (blue). H3K27me3 ChIP-Seq data was sourced from ENCODE. Negative control data (yellow) was generated using a Rabbit IgG antibody in CUTANA. Each dataset of CUTANA was acquired using 0.5 million cells.




FIGURE 3

Representative tracks of H3K4me3 on chromosome 7 from K-562 cells, comparing data generated from standard ChIP-seq methods and sequencing depth (black) to CUTANA (blue). H3K4me3 ChIP-seq data was sourced from ENCODE. Negative control data (yellow) was generated using a Rabbit IgG antibody in CUTANA. Each dataset of CUTANA was acquired using 0.5 million cells.

With CUTANA you can do more, faster, cheaper

CUTANA assays have low background and are highly reproducible! Cut through the noise with CUTANA




FIGURE 4

A representative 350 kb region of an H3K4me1 profile in K-562 cells, generated using CUTANA (yellow panels), native ChIP-seq (blue panels), or cross-linked ChIP-seq (green panels). All data were generated by EpiCypher and are expressed as reads per million (RPM). Color-coded gradient (to left) represents signal/noise (S/N) ratios determined by genome-wide analysis (bamFingerprint data, not shown).

Get inspiration to design your next experiment!

Check out these papers to see validated approaches for your ChIC/CUT&RUN workflows

ChIC workflows

- Schmid et. al, Mol. Cell 2004 (PMID : 15469830)
- Ku et. al, Nat. Methods 2019 (PMID : 30923384) ****Single Cell****

CUT&RUN workflows

- Skene and Henikoff, eLIFE 2017 (PMID : 28079019)
- Thakur and Henikoff, G&D 2018 (PMID : 29386331)
- Liu et. al, Cell 2018 (PMID : 29606353)
- Skene et. al, Nat. Protoc. 2018 (PMID : 29651053)
- Janssens et. al, Epi. Chromatin 2018 (PMID : 30577869)
- Brahma and Henikoff, Mol. Cell 2019 (PMID : 30554944)
- Oomen et. al, Genome Res. 2019 (PMID : 30655336)
- Zheng and Gehring, Plant Reprod. 2019 (PMID : 30719569)
- Ernst et. al, Nat. Commun. 2019 (PMID : 30890697)
- Hainer et. al, Cell, 2019 (PMID : 30955888)
- Meers et. al, eLIFE 2019 (PMID : 31232687)*
* Paper describes optimized protocol using pAG-MNase
- Meers et. al, Mol. Cell 2019 (In Press)

ORDERING INFO

CUTANA pAG-MNase for ChIC/CUT&RUN

Catalog No 15-1016

Coming Soon

CUTANA pAG-MNase assay kit
for ChIC/CUT&RUN

CUTANA pAG-Tn5
for ChIC/CUT&TAG

CUTANA pAG-Tn5 assay kit
for ChIC/CUT&TAG

Services

New to **CUT&RUN**? Lack the resources
to perform assay internally?

We can help!

Inquire at info@stratech.co.uk to learn
more about EpiCypher's 'end-to-end'
assay services for ChIC/CUT&RUN

www.stratech.co.uk/epicypher

+44 (0) 1638 782600

info@stratech.co.uk

