

FITC IgG Labeling Kit

For labeling 3 x 10 µg IgG

Reagent Storage

Kit Component	Storage Temp	Storage Notes
FITC Labeling Reagent	-20°C	Keep the vial in the desiccated container as supplied in the kit
Quenching Reagent B	-20°C or 2-8°C	Does not need to be kept desiccated.
100 mg/ml BSA (with 0.05% Sodium Azide)	-20°C or 2-8°C	Does not need to be kept desiccated.
Zeba Desalting Column with Collection Tube	2-8°C	Does not need to be kept desiccated. Do not store at freezing temperatures.

Introduction

The IgG FITC Labeling Kit utilizes a novel chemistry to generate highly reproducible fluorescein-labeled IgG with a simple procedure.

Features

- Pre-measured active labeling reagent.
- Can be used with up to 10 mg/ml (1%) BSA as a carrier protein.
- Completely scalable: conjugate anywhere from 0.1 to 1 gram IgG per reaction.
- Highly efficient FITC incorporation.
- Optimized FITC:IgG ratio by labeling time and temperature.
- Purification not usually required.

Products & Contents

Catalogue Number	FITC-Link-CC
For Labeling:	3 x 10 ug IgG
FITC Labeling Reagent	3 X 2 µL
Quenching Reagent B	100 µL
100 mg/ml BSA (with 0.05% Sodium Azide)	100 µL
Zeba Desalting Column with Collection Tube	3 each

Additional Reagents Required

None

Shelf Life

The performance of the product is guaranteed for a minimum of 12 months when stored as directed.

IgG Amount and Concentration and Buffers

The IgG to be labeled should be at a concentration 0.5 - 2.0 mg/ml in 1X PBS, pH 7.2 – 7.5. The IgG solution may contain up to 10 mg/ml (1%) BSA.

Conjugation Procedure for 0.1 mg of IgG

Note: Fluorescein is light sensitive. Perform reactions in the dark and keep light exposure to a minimum.

1. Let kit come to room temperature for at least 1 hour. Kit can be removed from -20°C up to 24 hours before use.
2. Spin the 1.5 ml screw-cap tube containing FITC Labeling Reagent for 1 minute at high speed ($> 14,000 \times g$).
3. Add 10 μ L IgG solution to the tube containing the FITC Labeling Reagent.
4. Vortex thoroughly for 10 seconds, then shake solution down to the bottom.
5. Incubate the labeling reaction at 37°C for 60 minutes in a heat block.

Note: Room temperature incubations will be successful. 15 minute – 4 hour incubation at 37°C will be successful.

6. Add 2 μ L of Quenching Reagent B and vortex.
7. Incubate at room temperature in dark for 5 minutes.
8. Desalt the complete reaction volume into pure 1X PBS using the included Zeba spin column. See the attached desalting protocol.
9. Store FITC-labeled IgG in the dark at 2-8°C.

Optional: Add 100 mg/ml BSA solution to achieve the desired final concentration of BSA.

Optional: Add glycerol to a final concentration of 40-50% (Not Provided).

Recommended Accessories

To remove excess FITC from the labeled IgG - Order from ThermoFisher:

Sample Size	Description	Cat #
2 – 12 µL	Zeba Spin Desalting Columns, Micro (75µL), 7K MWCO	89877, 89878
30 - 130 µL	Zeba Spin Desalting Columns, 0.5mL, 7K MWCO	89882, 89883
200 – 700 µL	Zeba Spin Desalting Columns, 2mL, 7K MWCO	89889, 89890
500 – 2000 µL	Zeba Spin Desalting Columns, 5mL, 7K MWCO	89891, 89892
700 – 4000 µL	Zeba Spin Desalting Columns, 10mL, 7K MWCO	89893, 89894

For concentrating IgG before labeling or for concentrating the final FITC-labeled IgG
– Order from MilliporeSigma:

Sample Size	Description	Cat #
Up to 500 µL	Amicon Ultra-0.5 Centrifugal Filter Unit with Ultracel-50 membrane	Z740176
Up to 2 mL	Amicon Ultra-2 Centrifugal Filter Unit with Ultracel-50 membrane	UFC205024
Up to 4 mL	Amicon Ultra-4 Centrifugal Filter Unit with Ultracel-50 membrane	UFC805008
Up to 15 mL	Amicon Ultra-15 Centrifugal Filter Unit with Ultracel-50 membrane	Z648000

Zeba Spin Desalting Column Instructions

75 μ L columns for 2 – 12 μ L Samples

Notes:

- Each column can desalt a 2-12 μ L sample.
- The resin slurry contains 0.03% sodium azide.
- These columns are recommended for desalting molecules > 7000 Daltons.

Storage:

- Store columns at 2-8°C.
- Columns may be stored at room temperature for several days without adverse effects.

Additional Materials Required

- Variable-speed bench-top microcentrifuge
- 1.5 mL microcentrifuge collection tubes

Spin Column Preparation

1. Remove column's bottom closure and loosen cap (do not remove cap).
2. Place column in the collection tube provided.
3. Centrifuge at 1000 $\times g$ for 1 minute to remove storage solution.
4. Discard storage solution from the collection tube. Remove column cap. It is not necessary to replace the cap for subsequent steps.
5. Add 50 μ L of buffer on top of the resin bed. Centrifuge at 1000 $\times g$ for 1 minute to remove buffer. Discard the buffer from the collection tube.
6. Repeat Steps 5 two additional times for a total of 3 column washes.

Sample Desalting

7. Place column in a new microcentrifuge tube and slowly apply the sample to the center of the compacted resin bed.
8. Do not add a stacker after application of the sample to the resin bed.
9. Centrifuge column at 1000 $\times g$ for 2 minutes to collect desalted sample. Discard desalting column after use.

Zeba Desalting Columns are a product of Thermo Fisher Scientific Inc.

Sample Size	Description	Cat #
2 – 12 μ L	Zeba Spin Desalting Columns, Micro (75 μ L), 7K MWCO	89877, 89878