



FuGENE 6

Quick Protocol

Preparing the FuGENE 6® Transfection Reagent

1. Before use, allow the vial of FuGENE® 6 Transfection Reagent to reach room temperature
2. Mix by inverting or vortexing briefly. If a precipitate is visible, briefly warm at 37 degrees C then cool to room temperature

General Transfection Protocol

1. To a sterile tube or U- or V-bottom plate add room temperature medium so that the final volume after adding the FuGENE 6® Reagent (in step 2) and DNA (in step 4) is 100µl.
2. For a 3:1 FuGENE® 6 Transfection Reagent to DNA ratio, add 6µl of FuGENE® 6 Transfection Reagent directly to medium, and mix immediately. For other ratios, consult table below.

Tube Label	Medium Final (ul)	FuGENE® 6 Reagent (ul)	DNA (ug)
2:1	100	4	2
3:1	100	6	2
4:1	100	8	2
6:1	100	12	2

3. Incubate the FuGENE® 6 / Medium mixture for 5 minutes at room temperature.
4. Add 2ug of DNA to the FuGENE® 6/Medium mixture (0.2-1.0 ug/ul) to a final volume of 100ul total. Vortex immediately,
5. Incubate complex at Room temperature for a minimum of 15 minutes. Up to 30 minutes
6. Add 2-10ul of the FuGENE® 6 Transfection Reagent/DNA mixture per well to a 96-well plate containing 100µl of cells in growth medium. Mix by pipetting or using a plate shaker. Return cells to the incubator for 24-72 hours.
7. Measure transfection efficiency using an assay appropriate for the reporter gene. For transient transfection, cells are typically assayed 24-72 hours after transfection.
8. See additional protocol information in Technical Manual available on www.fugene.com
9. For additional support please visit us at www.fugene.com

Learn more at: www.fugene.com