Other Required Materials

- Cryopreserved or freshly isolated T Cells, or CD3+ T Cells
- RPMI 1640 supplemented with 2 mM L-glutamine (or equivalent)
- Heat inactivated Fetal Bovine Serum (FBS)
- Penicillin/streptomycin
- Recombinant human IL-2 (rIL-2)
- 2-Mercaptoethanol
- Cell culture vessels
- Humidified CO₂ incubator or bioreactor

Part A) Activation & Expansion

1 Medium preparation

RPMI ● 10% FBS ● 1% Pen/Strep ● 20-50 μM 2-Mercaptoethanol ● 30-100 U/mL rIL-2
Or any other compatible culture media

2 HMP preparation

Resuspend Aim-Tconv HMP in the vial by vortexing.

3 HMP seeding

Calculate desired HMP seeding volume per well. Freshly prepare adequate HMP dilutions in culture medium for easy and accurate pipetting. Aliquot diluted HMP in to each well.

Recommendations

- ◇ T Cell seeding density : 2×10⁵ - 3×10⁵ cell per cm²
- ◇ HMP to PBMC ratio : 2:1
- ◇ HMP to T cell ratio : 1:1

Plate	Area (cm ²)	Cell/ well	Cell/ cm ²	1x HMP ratio	2x HMP ratio
96-well	0.32	8×10 ⁴	2.5×10 ⁵	2 μL	4 μL
48-well	0.95	2×10 ⁵	2.1×10 ⁵	5 μL	10 μL
24-well	1.9	4×10 ⁵	2.1×10 ⁵	10 μL	20 μL
6-well	9.5	2×10 ⁶	2.1×10 ⁵	50 μL	100 μL
T25	25	6×10 ⁶	2.4×10 ⁵	150 μL	300 μL
T75	75	2×10 ⁷	2.7×10 ⁵	500 μL	1000 μL

4 Cell seeding

Aliquot resuspended T cells. Gently mix T cells and HMP by pipetting up/down 3 to 5 times ensure HMP and cells are evenly distributed under microscope.

5 Co-culture

Incubate in a humidified 5% CO₂ incubator at 37°C. Monitor T cell morphology and confluency every other day. DO NOT disturb HMP-cell interaction in the first 4 days.

6 HMP cleanup

HMP will self-degrade by hydrolysis after around 14 days under culture environment. Alternatively HMP can be rapidly enzymatically degraded by adding the included HMP digesting buffer directly into the culture medium at 1X final concentration. After incubating at 37°C for 1-2 hours, the HMPs will fully degrade.

7 Restimulation

Restimulation is required for expansion longer than 10 days. T cell growth usually slows down after day 10 postinitial activation. Recommended HMP to cell ratio for restimulation is 0.5:1 to 1:1.

Component Descriptions

A suspension of hydrogel microparticles (HMP) made of chemically crosslinked dextran. HMPs sized 12-18 μm were coated with phospholipid bilayer, with membrane docked T cell activating signal panels.

- ◇ Supplied at 4 x 10⁷ / mL, in 1.0 mL PBS per vial
- ◇ A vial of HMP digesting buffer (10X, 1.0mL), containing enzymes to digest HMP.

Stability and Storage

- ◇ Shipped with ice, keep product refrigerated (2 to 8°C).
- ◇ Stable at 4°C for 9 months. Contents are sterile in unopened tube.
- ◇ Do not subject product to freezing, high temperature (greater than 40°C).

Part B) Cell Phenotype Characterization

1 Typical T cell morphology

- ◇ Resting T cells: smaller in size, round shaped
- ◇ Activated T cells: larger in size, irregular

2 Characterize T cells with FACS Accordingly:

	Marker	Type
Lineage	CD3+	Pan T
	CD4+	CD4 T
	CD8+	CD8 T
	CD45RA+ CD62L+	T _{naïve}
	CD45RA- CD62L+	T _{central memory}
	CD45RA+ CD62L-	T _{effector memory}
	CD45RA- CD62L-	T _{effector}
Activation	CD25	Activated T
	PD-1	Exhausted T

Things to note

When seeding T cells with HMP

Ensure T cells and HMP are evenly distributed to maximize interaction.

When you see T cell cluster

HMP and T cells tend to aggregate in the well center over time. Gently shake the culture plate to redistribute HMR and T cells.
Avoid unnecessary pipette mixing in the first 4 days, disturbing the HMP-cell clustering will cause suboptimal cell growth.

When to add new medium or split?

Monitor T cell growth periodically by performing T cell sampling and counting periodically. Supplement fresh medium or pass the cells to new culture vessels when


- ◇ Color turns orange yellow (acidic, ~ pH 6.5)
- ◇ T Cell grows to over 3 x 10⁶ cells /mL

When to add new medium or split?


Day 8-10 after previous stimulation.

Part C) Representative Data

Human Pan T cells co-cultured with Aim-Tconv at different time points

Representative bright-field micrographs showing the activation and expansion of pan T cells at various time points. Scale bar = 50 μm. Pan T cells were activated and expanded using Aim-Tconv HMP at 1:1 HMP:cell ratio. Red signals indicate the HMPs.

- 1 T Cell growth (n=4 T cells from healthy human donors) after activation by Aim-Tconv HMP and other popular T cell activation reagents. Expansion fold = viable cells at time t / initial number of PBMCs.
- 2 T Cell viability (n=4 T cells from healthy human donors) after activation by Aim-Tconv and other popular T cell activation reagents.
- 3 The expression of early stage marker (CD69) and late stage marker (CD25) indicating T cell activation. Unstimulated cell was used as gating.
- 4 T cell memory subtypes of PBMCs activated by Aim-Tconv and other popular T cell activation reagents.

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