NZYTech qPCR Selection Guide



NZYTech's qPCR products were developed for highly accurate real-time PCR compatible with a vast range of thermocyclers.

Please use the following selection guide to choose the product that is most appropriate for your equipment:

PASSIVE REFERENCE DYE

EQUIPMENT *	High	Low	None
Agilent			
AriaMX		8	
MX3000P™, MX3005P™, MX4000P™		8	
Applied Biosystems™	_		
	Φ.		
7000/7300/7700	8	•	
7500/7500 FAST 7900/7900HT/7900HT FAST	•	•	
QuantStudio™ 6, 7, 12k Flex/ViiA7™	•	•	
StepOne™/StepOne™ plus	8	•	
stepone /stepone plus	•		
Bio-Rad®			
CFX96™ /CFX384™			8
Opticon™/Opticon™ 2			8
Fluidigm [®]			
BioMark™		8	
Illumina®			
Eco™			8
Qiagen			
Rotor-Gene™ 3000			8
Rotor-Gene™ 6000			8
Rotor-Gene™ Q			8
Roche			
Lightcycler® 96			8
Lightcycler® 480			8
Lightcycler® Nano		7	8
(*) For other equipment please contact NZYTech (in	nfo@nzytech.com)		